
Validity and Reliability of the VO₂ Master Pro

Introduction

VO₂ analyzers are worn by athletes to measure performance metrics such as VO₂ max and metabolic equivalents. These metrics are used by coaches to design individualized training programs. VO₂ measurements can be obtained either directly by measuring the heat a subject generates within a small enclosed chamber (direct calorimetry), or indirectly through analysis of inspired and expired gases (indirect calorimetry). Indirect calorimetry is the preferred method for measuring athletic performance because it is simpler to implement and allows for maximal testing (VO₂ max) on a stationary bike or treadmill.

The advent of technology begets VO₂ analysis tools that go beyond the conventional lab test. Newer developments in technology have taken VO₂ analyzers off the lab bench and into the field. However, most portable units require large battery packs, usually housed in awkward backpacks or waistbands, that may impede the athlete during testing. A study by the University of IOWA found that increasing body weight with a backpack caused a detrimental effect on VO₂ [1]. At fixed treadmill speeds, subjects had an increase in VO₂ when loaded with a pack. The VO₂ Master Pro is a novel, patented VO₂ analyzer entirely self-contained to a small mask worn on the subject's face, eliminating the constraints of a backpack.

This study aimed to determine the validity and reliability of the VO₂ Master Pro when compared to a VacuMed metabolic simulator. The VacuMed simulator has been validated with the industry-standard Douglas bag method [2] and has been used to validate popular VO₂ analyzers such as the Cosmed K5 [3]. The simulator's design is based on concepts discussed by Huszczuk and Wasserman [4].

Twenty-four different set points of ventilation and metabolism were applied to two VM Pros, to examine accuracy and repeatability over a broad range. Bland-Altman plots, ordinary least products regression (OLP), mean percentage differences and typical percentage error were analyzed. A two-way random intraclass correlation (ICC) model was used to determine ICC coefficients to further determine reliability.

Method

Two VM Pros were used in three tests to determine intra and inter reliability. For each test, the device was connected to the flow conduit of a metabolic simulator (VacuMed Model 17056). Atmospheric pressure, temperature and humidity were recorded and used to determine the appropriate metabolic simulator mass flow rate of gas, as per manufacturer instructions [5]. Data was recorded at each test condition (table 1) after 30s of controlled input passed, to allow the metabolic simulator to reach stasis. Reported outputs were collected for a minimum of 10 breaths and averaged into a single recorded value.

Validity was examined using Bland-Altman plots, OLP regression, mean percentage differences and typical percent errors. Unit reliability was examined using intraclass correlation coefficients, minimum detectable change and typical error.

Table 1: Controlled Inputs of VacuMed Metabolic Simulator

Step	Tidal Volume (L)	Respiratory Frequency (Breaths Per Minute)	VO ₂ STPD (mL)
1	1.5	22	1000
2	1.5	30	1300
3	1.5	40	1600
4	1.5	50	2000
5	2.0	25	1500
6	2.0	35	1750
7	2.0	45	2100
8	2.0	55	3000
9	2.5	22	2000
10	2.5	30	2300
11	2.5	36	2600
12	2.5	40	3000
13	3.0	20	2500
14	3.0	30	3000
15	3.0	35	3600
16	3.0	40	4000
17	3.5	20	3000
18	3.5	30	3300
19	3.5	36	3600
20	3.5	40	4000
21	4.0	20	3500
22	4.0	25	3625
23	4.0	30	3750
24	4.0	35	4000

Results

Table 2: VO₂ Master Pro Intra – Inter Device Reliability Results

	Ventilation	VO ₂
ICC (Intra-Rater)	1.000 (0.995 - 1.000)	0.999 (0.996 - 1.000)
ICC (Inter-Rater)	1.000 (1.000 - 1.000)	0.999 (0.999 To 1.000)
Standard Deviation (%)	0.66	2.21
Typical Error Percent (%)	0.62	0.52
Minimum Detectable Change Intra-Rater (%)	1.72	2.34
Minimum Detectable Change Inter-Rater (%)	1.43	3.45

Intra and inter repeatability coefficients were greater than 0.999, demonstrating high reliability. Both ventilation and VO₂ had an excellent correlation between the metabolic simulator and the VO₂ Master Pro, with R² values of 0.9998 (p < 0.01) and 0.9972 (p < 0.01) respectively (figures 1 and 2). The confidence intervals of ventilation and VO₂ indicated no statistically significant bias, with their slopes passing through 1 and intercepts through 0. Mean percent differences were less than 1%.

Table 3: Validation Results

	R ²	Slope (95% CI)	Intercept (95 % CI)	Mean % diff (Min - Max)	Typical Percent Error
Ventilation	0.9998	0.9977 (0.9916 To 1.0039)	-0.2981 (-0.8668 To 0.2707)	-0.62 (-2.85 To 0.25)	0.466
VO₂	0.9972	0.9951 (0.9719 To 1.0183)	15.461 (-51.729 To 82.704)	0.22 (-5.05 To 3.91)	1.56

Bland-Altman plots demonstrated excellent agreeability and low deviation across the entire range, with a mean error of 0.5L/min [-1.31,0.332] for ventilation (figure 3) and 2mL/min [-94.0,98.15] for VO₂ (figure 4), given a 95% confidence interval.




Figure 1: Ventilation Reported by Simulator Vs VO2 Master Pro (R² = 0.9998, p < 0.01).




Figure 2: VO₂ Reported By Simulator versus VO₂ Master Pro ($R^2 = 0.9972$, $p < 0.01$).




Figure 3: Bland-Altman Plot of Ventilation with mean of 0.5L/min [-1.31, 0.332].




Figure 4: Bland–Altman Plot of VO₂ with a mean of 2mL/min [-94.0,98.15].

Discussion

This research is subject to several limitations. The VacuMed metabolic simulator is only capable of producing VO₂ inputs of up to four liters per minute. This study tested VO₂ inputs from one liter up to the maximum. The Vacumed 17056 has a series of fixed tidal volume selections from one to four liters, in half-liter increments. As such, this study was unable to test any tidal volume measurements between the half-liter steps, above four liters, or below one liter. In future studies, validation at greater tidal volumes, ventilations and metabolic rates should be evaluated. In addition, validation at resting metabolic rates and ventilations should be examined.

Conclusion

The VO₂ Master Pro demonstrated excellent agreement with the simulator throughout the entire tested metabolic range. The mean differences of 2mL/min for VO₂ and 0.5l/min for ventilation demonstrated no proportional or fixed bias, with confidence intervals crossing through zero and one respectively. Typical percent error was less than 2% for all metrics, easily achieving the 3% recommendation by the Australian sports institute [6]. Intra and inter reliability of the VO₂ Master Pro were excellent with ICC coefficients greater than 0.99, and minimum detectable changes less than 3.5% for VO₂ and 1.5% for ventilation. These results are comparable to popular VO₂ analyzers such as the Cosmed K5 [3].

References

- [1] M. L. Puthoff, B. J. Darter, D. H. Nielsen and H. J. Yack, "The Effect of Weighted Vest Walking on Metabolic Responses and Ground Reaction Forces," *Medicine & Science in Sports & Exercise*, pp. 746-752, 2006.
- [2] J. Bunn, J. L. Pittsley, S. V. Baker and J. Yates, "Assessment of accuracy of the Vacu-Med 17053 calibrator for ventilation, oxygen uptake, and carbon dioxide production," *Respiratory Care*, vol. 56, no. 4, pp. 472-476, 2011.
- [3] L. Guidetti, M. Meucci, F. Bolletta, G. P. Emerenziani, M. C. Gallotta and C. Baldari, "Validity, reliability and minimum detectable change of COSMED K5 portable gas exchange system in breath-by-breath mode," *PLOS ONE*, 2018.
- [4] A. Huszczuk, B. Whipp and K. Wasserman, "A respiratory gas exchange simulator for routine calibration in metabolic studies," *European Respiratory Journal*, pp. 465-468, April 1990.
- [5] VacuMed, Operating & Service Manual: Metabolic Calibrator / Simulator No. 17056, Ventura, CA: Author, 2017.
- [6] R. T. C. G. Australian Institute of Sport, *Physiological Tests for Elite Athletes*, 2nd ed., Champaign, United States: Human Kinetics Publishers, 2013.